

NEO PEDIATRIC COT

Quality technology for the benefit of the little ones



FEATURES

STRUCTURE

- Epoxy-coated steel tube frame.
- Gas piston height adjustment with a lever situated on the Undercarriage.

CONFORT

- Removable methacrylate cradle.
- Useful surface: 100%.

MANOEUVRABILITY

- Pneumatically-driven Trendelenburg/ reverse Trendelenburg position to $\pm 13^\circ$.
- Pneumatically-operated height adjustment, from 640 mm to 940 mm.
- Four 65 mm diameter swivel castors (two brake castors, two freewheeling castors) ensure

SAFETY

- Easy manoeuvrability of the cot.
- The NEO cot conforms to Directive 93/42/EEC.

Overall dimensions	92 x 57 x 102 cm
Height	64/94 cm.
Reverse trendelenburg	+ / - 13°
Cradle size	65 x 36 x 22-2- cm
Safe working load	15 kg.
Cot weight without accessories	21 kg.